


Location **South East, Surrey**
<https://www.freeadsz.co.uk/x-553804-z>

£99 for One OR £184 for Two - Save £thousands£

• The design of this hearing aid is inspired by nature. This hearing aid will particularly suit people with slight to moderate hearing loss. People who wear glasses and also have Behind the Ear hearing aids will especially benefit from an aid completely in the ear canal. People who do not wish to advertise that they are wearing hearing aids will like this aid, which cannot be seen unless looking straight into the ear. • Red for the right ear and blue for the left ear, means you can adjust and retain the settings for your right or left ear. Each hearing aid can be adjusted for volume sensitivity to suit exactly the requirements for each ear. Special price for two hearing aids • Or buy two for the Special Price and keep one as a spare. • Each hearing aid comes with 4 sizes of ergonomically designed silicone ear plugs that fit into the ear canal so that you can try different sizes and select the cap that is correct and most comfortable. The cap should be kept clean and a brush is provided for this purpose. The silicone cap should be replaced every few months and are available from us. • As the aid fits into your ear it has a nylon wire to make it easy to remove the aid. • This aid is a joy to wear – you will hardly notice it but will be able to hear what people are saying. Finally, you can again actively participate in conversations! The sound processor amplifies the volume of noise of any kind and is individually adjustable according to your needs. • The volume control is easy to adjust. Sound receiver volume is adjusted by removing the aid from your ear and using the screwdriver provided to gently reduce the volume to the minimum and then gradually increase the volume. Wear it for a day at the lower volume before turning the volume up very slightly and repeating the process. Try different settings until you find the most comfortable. When the volume is too high the aid will tend to whistle. • Note that the hearing aid must NOT be worn when swimming or showering as it must

NOT get wet. Also do not place the aid in an oven or expose the aid to excessive heat. Damp or extreme heat where with the melting of the wax, the extreme heat can cause the battery to short out and stop working. • We provide 3 months' supply of batteries at a very small cost. The very small cost is just the cost of the battery. The rest A10 is given and we supply more of the leading brands, Rayovac, Varta, Eveready, Duracell, etc. 4 pairs of 6 Rayovac Hearing Aid Batteries and a pair of 3 months' supply of your purchase request for the next 3 months' supply. • When you receive your hearing aid(s) if you change your mind about your purchase you will receive a full refund. If you return the goods in the original packaging undamaged within 30 days. If your hearing aid develops a fault, not caused by any mistreatment, within one year from purchase, please return it with proof of purchase and you will receive a brand new replacement at no cost to yourself. • Shipped free to all UK addresses. • My Story: I paid £1500 for an in-ear hearing aid a few years ago. After 2 years I accidentally dropped it into a hand-basin of water. It was out of warranty and unrepairable. A very expensive hearing aid! Then the NHS supplied



From UK- Completely In the
Ear Digital Hearing Aid 3
Months Batteries



From UK- Completely In the
Ear Digital Hearing Aid 3
Months Batteries



From UK- Completely In the
Ear Digital Hearing Aid 3
Months Batterless



From UK- Completely In the
Ear Digital Hearing Aid 3
Months Batteries



change your packaging in mistreatment a brand new 1500 for water. It w



Ear Digital Hearing Aid 3 Batterless



022
hps://www.seeed.cc/uk/x-5538



I'll refund you
 aid develops
 return it with pr
 ped free to ill
 r 2 years I acc
 ry expensive

Eel
 Digital
 Hearing Aid
 3

Batteries

<https://www.hearings.co.uk/x-553804-z>



in the good in
it, not caused
of purchase and
addresses. • My
ly dropped it in
ng and! Then the



Final
ny
u will receive
ny: I paid
hand-basin
HS supplied a
Months
Batteries

behind the ear hearing aid, but as I wear glasses I found that the frame and the hearing aid behind my ear constantly clashed and required adjustment. Also the tube from hearing aid to ear canal was uncomfortable on my ear and also very unsightly and my wife said it made me look older. Also, there was an echo. NOW, I wear a red in-ear mini in my right ear and the sound quality is excellent and I hardly know I am wearing a hearing aid.....and most people I speak with do not know either. BEST of ALL, my wife has stopped complaining that I have misunderstood what she has said or that she has to repeat herself or shout at me, because I am happy to wear the aid from morning to night. • HOW TO SAVE A FORTUNE: NHS hearing tests are free – if you think you have hearing loss in one or both ears ask your GP to arrange a hearing test. Then, you can accept the behind the ear hearing aids provided by the NHS or if you want in in-ear, less visible aid that will not clash with your glasses or sun-glasses purchase your new in-ear aids from us. The price at Boots, Specsavers, Hidden Hearing and other high-street providers could be £1000 or more. Technical Parameters: • The hearing aid is built to EU standards and bears the CE mark and also is certified that it contains no dangerous materials and has the ROHS mark. • MAX-SSPL90: 118+3 dBSPL • HFA-SSPL90: 110±4dBSPL • Peak Gain: 33dB • HFA Full-on Gain: 30±5dB • Frequency Range: 200-4000HZ • THD : ?10% • Eq. Input Noise Level: 28dBSPL • Battery Current Drain: 0.6mA Each Hearing Aid Package includes: 1 x hearing aid 4 x Different size earplugs 1 x screwdriver 1 x cleaning brush 1 user instruction leaflet 4 packs of 6 A10 Rayovac Hearing Aid