

New Blood Pressure monitor A amp D UA- 774 BPM Dual Memory (30 GBP)



Location **South West, Hampshire**
<https://www.freeadsz.co.uk/x-583268-z>



NEW A & D Medical UA-774 Blood Pressure Monitor has dual 30 measurement memories suitable for use by two people and has separate coloured illuminated start buttons.

Travel case included.

Two automatic one touch operation buttons.

Clinically validated.

Advanced oscillometric measurement.

World Health Organisation blood pressure classification indicator.

Last reading recall.

Large, easy-to-read, wide angle LCD 3-line display.

Displays systolic, diastolic and pulse readings simultaneously.

Fits bicep circumference 22-32cm.

Clock display with reminder alarm.

Automatic power-off.

Uses 4 AA batteries (included).

Standard AC adaptor inlet.

The UA-774 is clinically validated to the highest British Hypertension Society standard and uses an oscillometric method with professional accuracy. It is designed to detect the pulse and to automatically inflate the cuff to a systolic pressure level.

The Irregular Heartbeat Indicator lights up if an irregular heartbeat is detected during the measurement

and helps users to find an irregular heartbeat which may be associated with cardiovascular disease (this is for information only and is not a medical diagnosis). It is designed to detect the pulse and to automatically inflate the cuff to a systolic pressure level.

How to use:

Sit down comfortably and support the arm on a surface so the centre of the cuff is the same height as your heart. Relax for 5-10 minutes before taking the reading. Remain still and keep quiet during measurement. Do not measure right after physical exercise or a bath. Rest for 20-30 minutes before taking the measurement. During measurement, it is normal for the cuff to feel very tight. Refer to the included manual for descriptions of symbols displayed on the UA-774 once it has taken




New Blood Pressure monitor
A amp D UA- 774 BPM Dual
Memory

<https://www.freeadsz.co.uk/x-583268-z>



New Blood Pressure monitor
A amp D UA- 774 BPM Dual
Memory

<https://www.freeadsz.co.uk/x-583268-z>




New Blood Pressure monitor
A amp D UA- 774 BPM Dual
Memory

<https://www.freeadsz.co.uk/x-583268-z>



New Blood Pressure monitor
A amp D UA- 774 BPM Dual
Memory

<https://www.freeadsz.co.uk/x-583268-z>



New Blood Pressure monitor
A amp D UA- 774 BPM Dual
Memory

<https://www.freeadsz.co.uk/x-583268-z>



New Blood Pressure monitor
A amp D UA- 774 BPM Dual
Memory

<https://www.freeadsz.co.uk/x-583268-z>



New Blood Pressure monitor
A amp D UA- 774 BPM Dual
Memory

<https://www.freeadsz.co.uk/x-583268-z>



New Blood Pressure monitor
A amp D UA- 774 BPM Dual
Memory

<https://www.freeadsz.co.uk/x-583268-z>



New Blood Pressure monitor
A amp D UA- 774 BPM Dual
Memory

<https://www.freeadsz.co.uk/x-583268-z>



New Blood Pressure monitor
A amp D UA- 774 BPM Dual
Memory

<https://www.freeadsz.co.uk/x-583268-z>